Manufacturer

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ENGLISH

Composition

Geistlich Fibro-Gide® is a resorbable, porous collagen matrix obtained by standardized, controlled manufacturing processes. The collagen is obtained from veterinary certified pigs and is carefully purified to minimize the risk of immunological reactions. The device is weakly cross-linked to improve the volume stability of the material. Geistlich Fibro-Gide® is sterilized by irradiation. 

Major components of Geistlich Fibro-Gide® are Collagen Type I/III and Elastin.

Properties / Action

Geistlich Fibro-Gide® exhibits mechanical properties that enable soft-tissue augmentation. Weak cross-linking of the collagen strands is performed in order to improve the volume stability of the device and to allow the formation of new soft tissue. The device is fully degradable over time. In addition, Geistlich Fibro-Gide® provides sufficient space for the ingrowth of desired cells that generate new soft tissue and enables vascularization. The thickness as well as the other dimensions can be easily adapted to the desired size using surgical instruments such as a scalpel or scissors. Due to the good biological properties (good biocompatibility and low antigenicity) Geistlich Fibro-Gide® integrates well into the surrounding soft tissue. Due to the natural origin of the collagen a slight bulging and variations in the thickness of the Geistlich Fibro-Gide® matrix may be noted. Geistlich Fibro-Gide® is of a spongy consistency and readily absorbs fluids. Soaking Geistlich Fibro-Gide® in blood or saline will result in a small increase of the dimensions.

Intended Purpose

Geistlich Fibro-Gide® is intended to be used as an implantable device for regeneration and augmentation of soft tissue in oral and maxillofacial surgery.

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Indications

Geistlich Fibro-Gide® is indicated to be used in the surgical treatment of the following clinical conditions in the oral and maxillofacial area:

Defects and deficiencies of soft tissue associated with:

- > Presence of recessions*
- > Lack of soft-tissue volume

The defects and deficiencies may result from traumatic injuries, pathological conditions (congenital or acquired), medical treatments and therapies, or from history of personal lifestyle habits.

*Data on the use of the product in recession defects Miller Class III and IV is currently not available.

Patient Target Group

Geistlich Fibro-Gide® is intended to be used in patients who require surgical interventions aiming for regeneration and augmentation of soft tissue in the oral and maxillofacial area. Data on the use of the product during pregnancy or lactation, or in children are not available. Women who are pregnant or breastfeeding and children should not be treated with Geistlich Fibro-Gide®.

Intended User

Geistlich Fibro-Gide® is to be used by licensed dentists/physicians trained in surgical management of soft tissue in the oral and maxillofacial area.

Clinical Benefits

- > Geistlich Fibro-Gide® shows results comparable to connective tissue grafts in terms of esthetics and soft-tissue quantity and quality over time.
- > Using Geistlich Fibro-Gide® in recession coverage leads to results comparable to autologous tissue grafts including root coverage and reduction of recession depth and width.
- > When using Geistlich Fibro-Gide® as alternative to connective tissue graft (CTG), harvesting of autologous tissue from the palate can be avoided: An additional surgery site is avoided (including donor site morbidity and potential complications) and for the patient, the post-surgical pain and discomfort are lower.
- > Two thirds (66.7%) of patients prefer Geistlich Fibro-Gide® over CTG treatment.
- > Geistlich Fibro-Gide® can be applied for soft tissue augmentation simultaneously with guided bone regeneration (GBR) around implants to avoid a second surgery.
- > Clinical experience shows low incidence of wound healing complications with Geistlich Fibro-Gide®.

Instructions for use

The general principles of sterile handling and patient medication must be followed when using Geistlich Fibro-Gide®.

- > The defect is exposed and the usual surgical procedures in creating a properly prepared site are undertaken.
- > Geistlich Fibro-Gide® can be cut to the required size and shaped with sterile surgical instruments such as scalpels or surgical scissors. The trimming can be performed in either dry or wet state. It is recommended not to compress the dry matrix in order to not affect the po-

rous structure. Geistlich Fibro-Gide® can be hydrated for example with the patient`s blood or with saline. Immediately after wetting, the matrix will swell, which needs to be taken into account when cutting in dry state.

Dimensional increase due to swelling of the device: 3 – 12%

- > Geistlich Fibro-Gide® should be placed on the prepared wound bed.
- > Complete penetration of the matrix by blood and exudates enables the formation of a blood clot.
- > In case fixation of Geistlich Fibro-Gide® is needed, it can be sutured with standard suture techniques. Caution should be taken when suturing: It is recommended to place the sutures carefully with smooth pulling-forces and not too tight sutures in order to avoid ripping of the matrix.
- > Once the Geistlich Fibro-Gide® is in place it has to be completely covered in order to support the wound healing by primary intention.
- > Any tension on the soft tissue around/over Geistlich Fibro-Gide® should be avoided in order to avoid any wound dehiscence.
- > It is not recommended to use the matrix in more than one layer.
- > Geistlich Fibro-Gide® can be used in combination with dental implants and/or biomaterials (e.g. bone substitutes, barrier membranes) for guided bone regeneration.
- > The aftercare using mouth disinfectants containing chlorhexidine or containing iodine in combination with povidone or containing tea tree oil can be used. The product properties are not altered when used according to the Instructions for use. A discoloration might be observed. Local use of Geistlich Fibro-Gide® with further medicinal products, alcohol, disinfectants or antibiotics has not been studied.

Post-operative Care

The usual post-operative care and medication for soft-tissue augmentation procedures should be followed. Further prosthetic treatment should only be carried out after a healing period, to ensure complete tissue regeneration. In case of matrix exposure and bacterial contamination, rinsing with a bactericidal solution is standard practice.

In the rare event that matrix removal is necessary, the tissues adjacent to the matrix should be anaesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual matrix. Following careful reflection of the surrounding tissue, the remaining portion of the matrix can be excised, and the area debrided to remove any inflamed or infected tissue. In the post-operative phase, complete removal may not be possible since the product is intended to resorb over time.

Limitations of use

Contraindications

- > Geistlich Fibro-Gide® must not be used in the presence of an acute infection in the surgical area.
- > Geistlich Fibro-Gide® is derived from a porcine source and must not be used in patients with known sensitivity to porcine material or collagen allergies.

Precautions

Geistlich Fibro-Gide® should be used with special caution in patients with:

- > chronic infection at the surgical site
- > clinical conditions impairing tissue healing, especially:
 - a) uncontrolled metabolic diseases
 - b) prolonged corticosteroid therapy
 - c) immune impairment
- > medication and therapies impairing tissue healing, especially:
 - a) radiotherapy
- > heavy smoking

As these could lead to reduced tissue healing, one cannot expect the same rate of healing and overall success of the entire procedure as in otherwise healthy patients.

Geistlich Fibro-Gide® is a typical non-metallic diamagnetic material, in common with most human tissues. It cannot be heated or act as a magnet during the magnetic resonance (MR) examination. Geistlich Fibro-Gide® has not been specifically studied in the MR environment. Geistlich Fibro-Gide® may be used with metal hardware. In this case, special caution under the MR environment may be required.

The product is intended for use in a single patient in a single medical procedure.

Any product not required in the single procedure should be discarded since sterility and stability cannot be assured once the original primary packaging is opened.

Side effects

As the matrix is a collagen product, in very rare cases allergic reactions may occur.

Insufficient fixation of device might lead to displacement or loss of device.

Possible side effects that may occur with any surgery include swelling at the surgical site, flap sloughing, dehiscence, hematoma, bleeding, local inflammation/necrosis, and infection. These side effects could lead to reduced tissue healing, unplanned surgical intervention, loss of implant. Other surgery-related disturbances that may appear include tissue loss and pain.

Other information

Do not re-use

Do not re-sterilize

Use-by date

Temperature limit

Keep dry

Do not use if package is damaged

-  Keep away from sunlight
-  Caution
-  Consult instructions for use
-  Sterilized using irradiation
-  Batch code
-  Catalogue number
-  Recycling symbol
-  Medical Device
-  Authorized representative in the European Community / European Union
-  Importer
-  Sterile Barrier

Prior to use, the integrity of the sterile barrier system must be visually inspected for breaches.
Do not use the product if packaging is damaged or opened unintentionally.

Presentation

Geistlich Fibro-Gide® is packed in a sterile double packaging:

Geistlich Fibro-Gide® 1 matrix, size 15 mm x 20 mm x 3 mm

Geistlich Fibro-Gide® 1 matrix, size 20 mm x 40 mm x 3 mm

Geistlich Fibro-Gide® 1 matrix, size 15 mm x 20 mm x 6 mm

Geistlich Fibro-Gide® 1 matrix, size 20 mm x 40 mm x 6 mm

Summary of safety and clinical performance (SSCP)

A SSCP can be accessed on the European database on medical devices (EUDAMED) under the following URL*:

<https://ec.europa.eu/tools/eudamed>

*the SSCP can also be requested by sending an email to sscp@geistlich.ch

Patient Implant Card (PIC)

All information in regard of the provided PIC can be accessed under the following URL:

pic.geistlich-pharma.com

Notice: Any serious incident that has occurred in relation to the device should be reported to the Manufacturer (Geistlich Pharma AG) and the national competent authority.

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