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CE Mark

Composition

Geistlich Mucograft® / Geistlich Mucograft® Seal is a resorbable collagen matrix obtained by standardized, controlled manufacturing processes. The collagen is obtained from veterinary certified pigs and is carefully purified to minimize the risk of immunological reactions. Geistlich Mucograft® / Geistlich Mucograft® Seal is sterilized by γ -irradiation. 

Geistlich Mucograft® / Geistlich Mucograft® Seal has two structures, which form a bilayer product. The dense structure consists of compact collagen and has a smooth texture. A second structure consists of a thick, porous collagen spongy scaffold.

Major components of Geistlich Mucograft® / Geistlich Mucograft® Seal are Collagen Type I/III and Elastin.

Properties / Action

The low immunogenicity and proven biocompatibility favor the use of Geistlich Mucograft® / Geistlich Mucograft® Seal in oral surgeries. The fibrous microstructure of Geistlich Mucograft® / Geistlich Mucograft® Seal readily absorbs fluid. The collagen fibers of Geistlich Mucograft® / Geistlich Mucograft® Seal form a basic tissue matrix, which adheres well to the surrounding and underlying tissues. The specially designed porous structure of Geistlich Mucograft® is highly hydrophilic. This allows the immediate blood uptake to enable blood clot formation. The dense structure of compact collagen protects the wound in open healing situations and allows tissue adherence and favorable wound healing. The degradation of Geistlich Mucograft® / Geistlich Mucograft® Seal takes several months depending on the surgical site, vascularization and patient comorbidities. Due to the natural origin of the collagen a slight bulging and variations in the thickness of Geistlich Mucograft® / Geistlich Mucograft® Seal matrix may be noted. Fixation with sutures or pins is possible. The matrix is not radio-opaque.

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Intended Purpose

Geistlich Mucograft® / Geistlich Mucograft® Seal is intended to be used in closed and open healing situations to:

- > close wounds of oral mucosa
- > support wound healing and regeneration processes in case of oral mucosa defects and deficiencies.

The concerned oral mucosa defects and deficiencies may be surgically created or result from traumatic injuries, pathological conditions (congenital or acquired), medical treatments and therapies, or history of personal lifestyle habits.

Indications

Geistlich Mucograft® / Geistlich Mucograft® Seal is indicated to be used in the surgical treatment of the following clinical conditions:

- > defects and deficiencies of oral mucosa* associated with:
 - a) presence of recessions**
 - b) lack of keratinized tissue
 - c) scar tissue.
- > oral mucosa wounds* associated with:
 - a) bone regeneration procedures
 - b) lack of soft tissue.

*Clinical data is available for soft tissue defects up to size 40x30 mm.

**Data on the use of the product in recession defects Miller Class IV is currently not available.

Patient Target Group

Geistlich Mucograft® / Geistlich Mucograft® Seal is intended to be used in patients who require surgical interventions aiming for regeneration and augmentation of soft tissue in the oral and maxillofacial area. Data on the use of the product during pregnancy or lactation, or in patients below the age of 8 is not available. Women who are pregnant or breastfeeding should not be treated with Geistlich Mucograft® / Geistlich Mucograft® Seal.

Intended User

Geistlich Mucograft® / Geistlich Mucograft® Seal is to be used by licensed dentists/physicians trained in surgical management of soft tissue in the oral and maxillofacial area.

Clinical Benefits

- > In treatment of gingival recessions, Geistlich Mucograft® / Geistlich Mucograft® Seal leads to results comparable to autologous tissue considering recession coverage and recession depth. Thereby, donor site morbidity and potential complications of an additional surgery site are prevented.
- > Treatment with Geistlich Mucograft® / Geistlich Mucograft® Seal leads to effective and long-term stable gain of keratinized tissue around teeth and implants and in vestibuloplasty, comparable to autogenous tissue. Thereby, donor site morbidity and potential complications of an additional surgery site are prevented.
- > When using Geistlich Mucograft® / Geistlich Mucograft® Seal as alternative to autologous tissue, harvesting of autologous tissue can be avoided. Thereby, donor site morbidity and potential complications of an additional surgery site are prevented, and for the patient post-surgical pain and discomfort are lower.

- > Using Geistlich Mucograft® / Geistlich Mucograft® Seal leads to less scarring and better esthetic outcome than using free gingival graft including better texture and color match.
- > Patients report less pre-surgical anxiety and less post-surgical pain and discomfort when treated with GMG compared to treatment with autologous tissue. Five years after treatment, the patients still prefer the esthetics of sites treated with GMG over sites treated with autologous tissue.
- > Geistlich Mucograft® / Geistlich Mucograft® Seal allows for significantly shorter surgery time compared to treatment with autologous tissue.

Instructions for use

The general principles of sterile handling and patient medication must be followed when using Geistlich Mucograft® / Geistlich Mucograft® Seal:

- > The defect is exposed by means of (a) properly prepared flap(s) or tooth/teeth extraction and the usual surgical procedures in creating a properly prepared site are undertaken.
- > Geistlich Mucograft® / Geistlich Mucograft® Seal is cut to the required size and shape if necessary.
- > Geistlich Mucograft® / Geistlich Mucograft® Seal must remain dry during trimming and application to the defect. Geistlich Mucograft® Seal is designed to support mucosa wound healing after tooth extraction.
- > The Geistlich Mucograft® Seal spongy structure is identified by an indented stripe to assure the correct side is placed facing the bone. The compact structure must face outward, away from underlying bone. To achieve adhesion, the matrix is applied over the prepared site and held carefully in place for a period to allow penetration of the matrix by blood and exudates, which allows good adaptation and adhesion of the matrix to the underlying surface and facilitates the formation of a blood clot.
- > No compression of Geistlich Mucograft® / Geistlich Mucograft® Seal: The collagen matrix should remain uncompressed before, during and after surgery.
- > No tension around Geistlich Mucograft®: Any tension of the soft tissues around Geistlich Mucograft® should be avoided.
- > When using Geistlich Mucograft® Seal in bone regeneration procedures, the underlying bone deficiency must be filled with a space-maintaining material (e.g. Geistlich Bio-Oss® Collagen).
- > Geistlich Mucograft® / Geistlich Mucograft® Seal can be used both in closed and in open healing situations. When used for treatment of deficiencies of the oral mucosa associated with presence of recessions, Geistlich Mucograft® should be used completely submerged.
- > Fixation with sutures or pins is possible
- > The aftercare using mouth disinfectants containing chlorhexidine or containing iodine in combination with povidone can be used. The product properties are not altered when used according to the Instructions for use. A discoloration might be observed. Mouth disinfectants containing tea tree oil should be avoided. Local use of Geistlich Mucograft® / Geistlich Mucograft® Seal with further medicinal products, alcohol, disinfectants or antibiotics has not been studied.

Special instructions for use in periodontology

A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial infection as well as adequate oral hygiene. Therefore, prior to surgical intervention, patients must receive a hygiene phase - treatment, consisting of oral hygiene instructions, scaling and root planing, and occlusal adjustment when indicated.

Post-operative care

The usual post-operative care and medication should be given. Further prosthetic treatment should only be carried out after a healing period, to ensure complete soft-tissue regeneration. In the rare event that matrix removal is necessary, the tissues adjacent to the matrix should be anaesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual matrix. Following careful reflection of the surrounding tissue, the remaining portion of the matrix can be excised, and the area debrided to remove any inflamed or infected tissue. In the post-operative phase, complete removal may not be possible since the product is intended to resorb over time.

Limitations of use

Contraindications

- > Geistlich Mucograft® / Geistlich Mucograft® Seal should not be used in the presence of an acute infection in the surgical area.
- > Geistlich Mucograft® / Geistlich Mucograft® Seal should not be used in patients with a known allergy to collagen.

Precautions

Geistlich Mucograft® / Geistlich Mucograft® Seal should be used with special caution in patients with:

- > chronic infection at the surgical site
- > clinical conditions impairing tissue healing, especially:
 - a) uncontrolled metabolic diseases
 - b) prolonged corticosteroid therapy
 - c) immune impairment
- > medication and therapies impairing tissue healing, especially:
 - a) radiotherapy
- > heavy smoking

As these could lead to reduced tissue healing, one cannot expect the same rate of healing and overall success of the entire procedure as in otherwise healthy patients.

Geistlich Mucograft® / Geistlich Mucograft® Seal is a typical non-metallic diamagnetic material, in common with most human tissues. It cannot be heated or act as a magnet during the magnetic resonance (MR) examination. Geistlich Mucograft® / Geistlich Mucograft® Seal has not been specifically studied in the MR environment.

Geistlich Mucograft® / Geistlich Mucograft® Seal may be used with metal hardware. In this case, special caution under the MR environment may be required.

Intraoperatively, if there is need to remove the product, complete removal can be achieved. In the post-operative phase, complete removal may not be possible since the product is integrated over time with the existing tissue.

The grafted area should not be brushed for the first four weeks and excessive muscle tractioning or trauma to the treated areas should be avoided for the first four weeks.

The product is intended for use in a single patient in a single medical procedure. Any product not required in the single procedure should be discarded since sterility and stability cannot be assured once the original primary packaging is opened.

Side effects

As the matrix is a collagen product, in very rare cases allergic reactions may occur. Insufficient fixation of device might lead to displacement or loss of device. Possible side effects that may occur with any surgery include swelling at the surgical site, flap sloughing, dehiscence, hematoma, bleeding, local tissue inflammation, local tissue necrosis, pain and infection. Incompatibility reaction with locally and systemically applied antibiotics cannot be excluded.

These side effects could lead to reduced tissue healing, healing disturbances, unplanned surgical intervention, loss of implant.

Other information

Do not re-use

Do not resterilize

Use-by date

Temperature limit

Keep dry

Do not use if package is damaged

Keep away from sunlight

Caution

Consult instructions for use

Sterilized using irradiation

Batch code

Catalogue number

Recycling symbol

Medical Device

Authorized representative in the European Community / European Union

Date of manufacture

Importer

Double sterile barrier system

Prior to use, the integrity of the sterile barrier system must be visually inspected for breaches. Do not use the product if packaging is damaged or opened unintentionally.

Presentation

Geistlich Mucograft® / Geistlich Mucograft® Seal is packed in a sterile double packaging:

Geistlich Mucograft® 1 matrix, size 15 mm x 20 mm

Geistlich Mucograft® 1 matrix, size 20 mm x 30 mm

Geistlich Mucograft® 1 matrix, size 30 mm x 40 mm

Geistlich Mucograft® Seal, 1 matrix, 8 mm diameter

Geistlich Mucograft® Seal, 1 matrix, 12 mm diameter

Summary of safety and clinical performance (SSCP)

A SSCP can be accessed on the European database on medical devices (EUDAMED) under the following URL:

<https://ec.europa.eu/tools/eudamed>

*the SSCP can also be requested by sending an email to sscp@geistlich.ch

Patient Implant Card (PIC)

All information in regard of the provided PIC can be accessed under the following URL:

pic.geistlich-pharma.com

Notice: Any serious incident that has occurred in relation to the device should be reported to the Manufacturer (Geistlich Pharma AG) and the national competent authority.

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